

DESCRIPTION

Densely striated, extremely strong tendons responsible for flexion of the foot

BRAND

FlexiGRAFT[®]

CLINICAL APPLICATIONS

- Anterior Cruciate Ligament Reconstruction
- Posterior Cruciate Ligament Reconstruction
- Tendon soft tissue augmentation



FEATURES AND BENEFITS

- Physician can predetermine diameter size of tunnels
- Physician can predetermine graft length for procedure
- Allows for various fixation techniques

ALTERNATIVE PRODUCTS

Measuring Tendons

The diameter (mm) is determined by pulling the folded tendon through the appropriate hole in the sizing block

Frozen	Description	Size
FPLT	Peroneus Longus Tendon	L ≥ 23 cm ; D ≥ 7.5 mm
FPLT-SL	Peroneus Longus Tendon: Short Length	L = 170-225 mm ; D ≥ 7.5mm