

DESCRIPTION

Part of the group of hamstrings/fleshy muscles and associated tendons located in posterior area of thigh

BRAND

FlexiGRAFT[®]



CLINICAL APPLICATIONS

- Anterior Cruciate Ligament Reconstruction
- Posterior Cruciate Ligament Reconstruction
- Tendon Soft Tissue Augmentations

FEATURES AND BENEFITS

- Physician can predetermine diameter size of tunnels
- Physician can predetermine graft length for procedure
- Allows for various fixation techniques

ALTERNATIVE PRODUCTS

Measuring Tendons

The diameter (mm) is determined by pulling the folded tendon through the appropriate hole in the sizing block

Frozen		Description	Size
FROPE		Semi Tendinosis or Gracilis	Length > 150 mm and Diameter = 4 to 5.5 mm
FST		Semi-tendinosus Tendon (Length = 160 - 250 mm)	L >= 230 mm D >= 4 mm
FSTP		Semi-tendinosus Tendon	L = 160 - 180 mm D = 4 - 6 mm