

DESCRIPTION

Cancellous or trabecular bone that is crushed - particle size varies - primarily from distal femur, condyles, proximal tibia, femoral head, distal tibia and humeral head.

BRAND

ReadiGraft® BLX

CLINICAL APPLICATIONS

- Fracture Management
- Filling Osseous Defects
- To Augment Total Joint Arthroplasty



FEATURES AND BENEFITS

- Provides an osteoconductive matrix for incorporation
- Allows physician to supplement with additional graft material
- Allows physician to use for various clinical applications

ALTERNATIVE PRODUCTS

Corticocancellous, Femoral Head, Humeral Head for grinding

Freeze-dried	Description	Size
DCAN5	Demineralized Cancellous, 1 - 8 mm	5 cc
DCAN40	Demineralized Cancellous, 1 - 8 mm	40 cc