

DESCRIPTION

Parallel cross-sections processed from either the humerus or tibia

BRAND

VertiGraft[®]

CLINICAL APPLICATIONS

- Anterior Lumbar Interbody
- Anterior Thoracic Fusion
- Corpectomy



FEATURES AND BENEFITS

- Cortical bone provides biomechanical strength
- Provided in a variety of sizes to accommodate surgeon's needs and preferences

ALTERNATIVE PRODUCTS

Iliac Crest Wedge, Patellar Wedge

Freeze-dried	Description	Size
VHUMCS8	Humeral Spacer Parallel	H = 8 mm D = 18 - 22 mm
VHUMCS10	Humeral Spacer Parallel	H = 10 mm D = 18 - 22 mm
VHUMCS12	Humeral Spacer Parallel	H = 12 mm D = 18 - 22 mm
VHUMCS14	Humeral Spacer Parallel	H = 14 mm D = 18 - 22 mm
VHUMCS16	Humeral Spacer Parallel	H = 16 mm D = 18 - 22 mm
VHUMW11	Humeral Spacer Lordotic	Ha = 14 mm Hp = 11 mm D = 18 - 22 mm
VHUMW13	Humeral Spacer Lordotic	Ha = 16 mm Hp = 13 mm D = 18 - 22 mm
VHUMW15	Humeral Spacer Lordotic	Ha = 18 mm Hp = 15 mm D = 18 - 22 mm
VHUMW17	Humeral Spacer Lordotic	Ha = 20 mm Hp = 17 mm D = 18 - 22 mm