



1864 Concert Drive, Virginia Beach, VA 23453

888.847.7831 | [www.AccessLifeNetHealth.org](http://www.AccessLifeNetHealth.org)

## DESCRIPTION

Tricortical graft processed from the calcaneus, also known as POD (Podiatry) Wedge



## BRAND

MatriGRAFT®

## CLINICAL APPLICATIONS

- Ankle Fusion Osteotomy
- Subtalar Joint Fusion
- Anterior Cervical Fusion

## FEATURES AND BENEFITS

- Allows for some structural support
- Allows physician to select appropriate graft for clinical procedure
- Provides an osteoconductive matrix for incorporation
- Pre-sized for convenience

## ALTERNATIVE PRODUCTS

Tricortical Patellar Wedge, Iliac Crest Wedge

Freeze-dried	Description	Size
POD.WEDGE	Calcaneus Wedge	H = 15 mm L >= 25 mm
POD.WEDGE8	Calcaneus Wedge	H = 18 mm L >= 25 mm