


DESCRIPTION

Cloward Dowel

Cylindrical graft with cortical face and cancellous body with pinset on top; processed from distal femoral condyle, proximal tibia, distal tibia or femoral head



BRAND

MatriGRAFT[®]

CLINICAL APPLICATIONS

- Anterior Cervical and Lumbar Fusion

FEATURES AND BENEFITS

- Allows for some structural support
- Allows for impactor to be screwed into graft
- Allows physician to select appropriate graft for clinical procedure
- Provides an osteoconductive matrix for incorporation
- Pre-sized for convenience

Freeze-dried	Description	Size
CL 10	Cloward Dowel	D = 10 mm L = 15 mm min
CL 11 (M18HS)	Cloward Dowel	D = 11 mm L = 15 mm min
CL 12	Cloward Dowel	D = 12 mm L = 15 mm min
CL 13 (M17HS)	Cloward Dowel	D = 13 mm L = 15 mm min
CL 14	Cloward Dowel	D = 14 mm L = 15 mm min
CL 16	Cloward Dowel	D = 16 mm L = 15 mm min
CL 18	Cloward Dowel	D = 18 mm L = 15 mm min
MR 11	Cloward Dowel	D = 11 mm L = 15 mm min; no septil hole

Preservon	Description	Size
-----------	-------------	------

PCL12

Cloward Dowel

Diameter = 12 mm